



## 2X PCR Bestaq™ MasterMix

Store at -20°C

Cat. No.	Description	Quantity
G464	2X PCR Bestaq™ MasterMix	5 ml
G464-dye	2X PCR Bestaq™ MasterMix with dye	5 ml

### Product Description

**crb's** 2X PCR Bestaq™ MasterMix is a ready-to-use mixture containing **crb's** Bestaq™ DNA Polymerase, dNTPs, and reaction buffer with proprietary additives in a 2X concentration. It contains all the necessary reagents for an efficient, high-yield amplification of DNA. The 2X PCR Bestaq™ MasterMix with dye (Cat. No. G464-dye) contains a green dye blend which resolves during gel electrophoresis into a turquoise band at ~4000 bp and a yellow band at the ~50 bp region. The resolving bands of dye allow easy visualization of the real time progress of your gel electrophoresis, and the yellow band indicates the migrating front in the gel.

**crb's** Bestaq™, with its superior fidelity and extreme robustness, is one of the best-performing DNA polymerases available on the market. With its monomeric, unique topological structure, Bestaq™ DNA Polymerase possesses enhanced processivity, elevated yield, boosted speed, extended amplification length, and improved difficult template tolerance. This high-performing polymerase and its MasterMix are the ideal choice for any PCR application, especially for cloning difficult templates, amplifying A/T- and G/C-rich sequences, and replicating long amplicons. With over 50X more accuracy than Taq DNA polymerase, Bestaq™ DNA Polymerase and its MasterMix can consolidate all PCR protocols and reactions into one simple, efficient, and extremely reliable PCR system. Bestaq™ DNA Polymerase generates blunt-end PCR products and has both 5'-3' polymerase activity and 3'-5' proofreading exonuclease activity.

### Features and Benefits

- Saves preparation time by combining Bestaq™ DNA Polymerase, dNTPs and reaction buffer into a ready-to-use mixture.
- Reduces the risk of contamination by decreasing the number of pipetting steps.
- Provides consistent reaction performance and results.

### Shipping and Storage

Keep at -20°C for long-term storage. 2X PCR Bestaq™ MasterMix and 2X PCR Bestaq™ MasterMix with dye are stable at 4°C for one month or for fifteen freeze-thaw cycles. For daily use, it is recommended to keep an aliquot at 4°C.

### Protocol

All PCR experiments should be assembled in a nuclease-free environment. In addition, DNA sample preparation, reaction set-up and subsequent reaction(s) should be performed in separate areas to avoid cross-contamination. The use of "clean" pipettors designated for PCR and aerosol-resistant barrier tips are recommended. Always keep the control DNA and other templates to be amplified isolated from the other components.

A negative control reaction (omitting template DNA) should always be performed in tandem with sample PCR to confirm the absence of DNA contamination.

1. Add the following components to a sterile 0.2 ml PCR tube sitting on ice.

Components	Volume	Final Concentration
Template DNA	~100 ng	~2 ng/μl
Forward primer (10 μM)	1 - 2.5 μl	200 - 500 nM
Reverse primer (10 μM)	1 - 2.5 μl	200 - 500 nM
2X PCR Bestaq™ MasterMix/ with dye	25 μl	1X
Nuclease-free H <sub>2</sub> O	up to 50 μl	-

- *We recommend preparing a mastermix for multiple reactions to minimize reagent loss and enable accurate pipetting.*
2. Mix contents of tube and centrifuge briefly.
  3. Incubate tube in a thermal cycler at 94°C for 3 mins to completely denature the template.
  4. Perform 30 - 35 cycles of PCR amplification as follows:
    - Denature:** 94°C for 10 secs
    - Anneal:** 45 - 72°C for 30 secs
    - Extend:** 72°C for 1 min/3 - 4 kb template
  5. Incubate for an additional 5 mins at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.
  6. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide or SafeView™ (Cat No. G108) staining. If 2X PCR Bestaq™ MasterMix with dye is used, load the samples directly without adding additional loading dye. Use appropriate molecular weight standards.



All **crb** PCR, RT-PCR, and qPCR products are ISO 13485:2003 and 13485:2012 certified as diagnostic grade and in compliance with all regulatory requirements for the design and manufacture of medical devices, as outlined by the International Organization for Standardization (ISO). For technical questions, please email us at [support@coderegenesis.com](mailto:support@coderegenesis.com) or visit our website at [www.coderegenesis.com](http://www.coderegenesis.com).